General remarks

Cholesterol in human body is present as part of cellular membranes, is pre-cursor of steroid hormones and of bile acids; for the most part it comes from liver and from feeding. Cholesterol is transported in plasma by lipoproteins, complexes formed by lipids and apo-lipoproteins. Prevention and early diagnosis of the cardiovascular diseases are mainly based on determination of hypercholesterolemia as a risk factor.

Using Visual Cholesterol test it is possible to check quickly cholesterol level in blood, so to contact a physician in case that level is high and in this way to significantly reduce heart disease risk.

Cholesterol level in blood is affected by diet and also by: drugs, diabetes mellitus, pregnancy and serious diseases. Hypercholesterolemia could be a secondary effect of other disease (thyroid insufficiency, biliary cirrhosis, pancreatic disease), or essential. Perform the test not before the physiologic conditions are normal (for ex. 3 months after childbirth). Perform test at fast from at least 6-9 hours.

Procedure

Prepare a chronometer or a clock with seconds. Take the comparative colour card outside the pack.

Open vial and take the reactive strip soon before executing the test. Vial and cap can be eliminated together with normal domestic waste.




Fingerpricking

Read instructions carefully.

Wash hands with hot water and soap and dry well them. If disinfectants are used, before performing puncture of fingertip, be sure that this part is dried.

Follow pictures:

- Be sure that protective cap is not removed. Load lancet device pushing protective cap inside it.
- Twist off the protective cap as described in the second picture and then take out when it easily turn.
- Press the disinfected fingertip on the part of lancet device where protective cap is just removed, be careful to position it not too close the nail and not in the centre of fingertip. Then press little to liver to hear a click. In the case the puncture was not performed, it is possible to re-load the lancet device housing-back the needle by pressing the protective cap into his original site. Lancing device is single use device: do not use to perform a second puncture.
- Let go out drop of blood and put it on the white reactive area in the centre of the green zone without touching the reactive area with the finger.
- Wait about 35-40 seconds to let the chemical reaction to develop.
- Turn the strip on the backside.
- Within next 25 seconds compare the developed colour with each colour area of the comparative card by positioning of the backside of the strip behind the comparative card in correspondence with the holes. Read on the comparative card the written respective cholesterol concentration.

Interpretation of results

Since the test result is a semi-quantitative evaluation involving colour, shade and intensity, it is useful the opinion of a second person on final result. It is important to perform the colour comparison within next 25 seconds, otherwise the colour will going to be lighter and the evaluated concentration will be lower than the real one. Each colour area corresponds to a cholesterol value. Lightest colour is cholesterol 130 mg/dl. A reactive strip with a colour lighter than 130 mg/dl has to be considerate as a cholesterol level lower than 130 mg/dl. A reactive strip with a colour darker than 350+ mg/dl has to be considerate as a cholesterol level higher than 350 mg/dl.

Normal and abnormal values (according to the latest criteria)

- Less than < 200 mg/dL = Desirable Blood Cholesterol
- 200 to 239 mg/dL = Borderline High Blood Cholesterol
- Greater than > 240 mg/dL and above = High Blood Cholesterol

Important! All results higher than 200 mg/dl have to be repeated with a laboratory quantitative method before to make a hypercholesterolemia diagnosis.

Principle of the test

Test is based on the reaction Cholesterol Esterase/Cholesterol Oxidase/Peroxidase/chromogen. Intensity of the colour developed by the reaction is proportional to cholesterol concentration in the blood.

Limitation

- * System is optimised for use of capillary blood. Do not use serum or plasma.
- * Reducing substances as ascorbic acid in normal concentration do not affect results, but if assumed in high quantity could give final colour lighter than real one.
- * High haemoglobin and/or bilirubin levels could interfere with result determining a lighter colour than the real one.

* Colour-blind person shall not use comparative colour procedure.

* High hematocrite value (>50%) gives a lighter colour and a cholesterol concentration lower than the real one, while a low hematocrite value (<30%) gives a darker colour and a cholesterol concentration higher than the real one.

* Result is affect by high levels of steroids.

Warning

Strips have to be stored between 10°C and 35°C (50°F e 95°F) and in their original capped vial. Do not expose to direct sunlight. Do not use within and not after expiry date indicated on the pack and within 3 months from first opening.

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