Strips for semi-quantitative determination of glucose in whole blood

General remarks

The patient can perform the test by himself; this must however be done under the constant control of a physician or a medical staff to define together the timetable of the determinations, the interpretation of results and the therapy.

Description

Visual Glucose permits the determination of glucose concentration in human capillary blood. Glucose concentration is estimated by comparison with comparative card in the range 25-500 mg/dl. Comparative card is attached in the pack.

Indication

For rapid semi-quantitative determination of blood glucose in diabetic person or in diagnosing suspected diabetics, in suspected hypoglycaemic states and for diabetic self-monitoring.

Principle

Each strip has a reactive area on which Glucose-oxidase reacts with D-glucose present in blood and generates hydrogen peroxide which, due to the presence of peroxidase, reacts with the chromogen. The colour intensity is proportional to the concentration of glucose present in the blood drop.

Test Preparation

Prepare a chronometer or a clock to read seconds. Pull out comparative card from pack. Open the vial and pull out reactive strip immediately before execution of the test. Vial and cap can be eliminated together with normal domestic waste.


Executing the test

Required accessories, not included: lancet device and sterile lancets.


Finger pricking

Read instructions carefully

1. Wash hands with hot water and soap and dry well them.
2. Perform puncture of fingertip
3. Remove the first blood drop with sterile gauze
4. Press the finger tip to obtain a second blood drop and approach the finger to strip. When there is a hanging drop (about 10µl), put it on the white reactive area in the centre of the blue zone without touching the reactive area with the finger.
5. Be sure that blood drop completely cover reactive area of the strip
6. Wait about 50 seconds to let the chemical reaction to proceed.
7. Turn the strip on the backside
8. Within next 25 seconds compare the developed colour on the strip with each coloured area of the comparative card by positioning the backside of the strip on the rear of the comparative card in holes. Find the closest matching colours. Read on the comparative card the written respective glucose concentration.

50"

25"

Visual reading is affected by personal evaluation and can change individual by individual.

Interpretation of results

Since the test result is a semi-quantitative judgement based on colour, shade and intensity, it is useful to get the opinion of a second person on final result. It is important to make the colour comparison within next 25 seconds, otherwise the colour will going to be darker and the evaluated concentration will be more than the real one. Each colour block corresponds to a glucose value. Lightest colour is glucose 25 mg/dl. A reactive strip with a colour lighter than 25 mg/dl has to be considered as a glucose level lower than a 25 mg/dl. A reactive strip with a colour darker than 500+ mg/dl has to be considered as a glucose level higher than a 500 mg/dl.

Normal blood glucose values

At fasting normal glycaemia values in a non-diabetic person are: 60-100 mg/dl (3,3 - 5,6 mmol/l) on whole blood or 70-110 mg/dl (3,9 - 6,1 mmol/l) on serum or plasma.

Important! Caution is recommended in the interpretation of glucose values below 50 mg/dl or above 250 mg/dl. Consult a physician as soon as possible if such values are obtained.

Limitation

- * The system is intended for analysis with capillary blood. Do not use serum or plasma.
- * If a glucose value is obtained which does not correspond to the sensations felt by the patient and to the knowledge he usually has of his health condition and symptoms, repeat the test. Eventually, consult your physician.
- * Capillary blood glucose values may be significantly lower than "true glucose levels" in the hyperglycaemic hyperosmolar state, with or without ketosis. Critical patients should not be tested by this system or at least tested with extreme caution.
- * Hands and fingers contaminated with sugar from foods or beverages may cause false elevated results. Wash hands in warm soapy water, rinse and dry before testing.
- * Severe dehydration and excessive water loss may cause inaccurately low results. Dehydrating factors include vomiting, diarrhoea and the intake of certain drugs such as diuretics.
- * Reducing agents e.g. ascorbic acid, at normal concentrations, do not interfere. If high quantities are intaken, abnormally low results could be obtained.
- * Grossly lipemic samples may interfere.
- * Colour-blind individuals should not use visual colour procedures.
- * Abnormal hematocrit values (higher than 55% or lower than 35%) can affect the results with visual glucose strips.
- * Visual glucose strips are not suitable to perform determinations on neonatal blood samples.

Warning

Strips have to be stored between 10°C and 35°C (50°F and 95°F) and in their original capped vial. Do not expose to direct sunlight. Do not use within and not after expiry date indicated on the pack and within 3 months from first opening.

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